libShower

Manual V0.1 frank.meisel@lhep.unibe.ch

Description

 This manual is intended to make the user familiar on the way showers are dealt with

The user should not bother about details of programming


 To increase the performance, the user may give additional information on the date to obtain a better result

Overview

- Shower for Event classification.
- Shower characteristics
- Shower implementation in FEDRA


Object description

- A track in fedra
 - One Segment Per Plate
 - N(), Npl(), N0, Theta(), XYZ, P()



Object description

- A shower in fedra
 - Can have More than One Segment Per Plate
 - N(), Npl(), N0, Theta(), XYZ, P()
 - Plus additional information: Iongitudinal&transversal Profile!



Object description

- Using EdbTrackP as object class for Shower is ok, but we lose physical information.
- Compensated by infos in EdbShowerAlgE...
 - May be put later in a extra EdbShowerP class...
- Up to now: how is a shower represented/stored?
- We have Shower.root which uses a TTree to store the variables. Technically equivalent to do storage like EdbTrackP, but it requires complex I/O conversions to get from
 - TTree entry → EdbTrackP

What is to be done for showers?

- In case we have scanback done, vertex found, scanforth done (hadr. re-interaction), we need for the full description of the event some other information: showers?
- Why are showers for the event classification important?
- Electrons: nu_e from beam; tau->e;
- Photons: present in almost any event, mainly through the decay of:
- Pi0: Find two photons, match them in case of more

Finding showers

- Showers are either attached to vertex directly or indirectly:
 - Electrons: start to shower direct from the vertex without flight length; IP to vertex rather small; DeltaZ to vertex very small.
 - Photons: fly then showering happens; therefore IP to vertex bigger than attached tracks; DeltaZ to vertex can be large (mean flight length ca10plates).
- Finding a electron shower can be easier than photon shower (directly attached track or BT)

Properties of Showers

- We need to be sure of the following things.
- We have to find the shower (start)
- We have to reconstruct (collect all the tracks) it having as much of it reconstructed as possibles (lowest loss as possible).
- After shower(s) found and reconstructed we need the main properties:
 - Energy
 - ID
 - 1ry-2ry vtx attachment

Properties of the Shower.root file. a) basic quantities

Name	in	Shower		root
------	----	--------	--	------

Description

Analogy in EdbTrackP

->MCEvt()

number eventb/l

sizeb nfilmb[] lengthfilmb

-b[] X, V, Ztx,ty -b[]

plateb -b[] MC event number

Number of Basetracks

Plate ID of BT[i]

Number of Plates Crossed

X,Y,Z Coordinates of BT[i]

TX,TY Coordinates of BT[i]

??(to be looked up)

N()

??(to be looked up)

Npl()

->X(),Y(),Z()

->TX(),TY()

??(to be looked up)

Numbereventb[]

showerID/I

isizeb/l

ntrace1simu[sizeb]/I

ntrace2simu[sizeb]/I

ntrace3simu[sizeb]/I

Ntrace4simu[sizeb]/I

chi2btkb[sizeb]/F

Numbereventb[]

showerID/I

isizeb/l

ntrace1simu[sizeb]/I

ntrace2simu[sizeb]/I

ntrace3simu[sizeb]/I

Ntrace4simu[sizeb]/I

chi2btkb[sizeb]/F frank.meisel@lhep.unibe.ch

Numbereventb[]

showerID/I

isizeb/l

ntrace1simu[sizeb]/I

ntrace2simu[sizeb]/I

ntrace3simu[sizeb]/I

Ntrace4simu[sizeb]/I

chi2btkb[sizeb]/F

10

Properties of the Shower.root file. b) derived quantities

Name in Shower . root

Description

Analogy in EdbTrackP

number eventb/l

number eventb/l

number eventb/l

sizeb nfilmb[] lengthfilmb -b[] X, V, Z

sizeb nfilmb[] lengthfilmb

x,y,z-b[]

tx,ty -b[] plateb -b∏

lengthfilmb X, Y, Z-b[]

sizeb

nfilmb[]

tx,ty -b[]

plateb -b[]

plateb -b∏

Numbereventb[]

showerID/I

isizeb/l

tx,ty

ntrace1simu[sizeb]/I

-b[]

ntrace2simu[sizeb]/I

ntrace3simu[sizeb]/I

Ntrace4simu[sizeb]/I

chi2btkb[sizeb]/F

Numbereventb[]

showerID/I

isizeb/l

ntrace1simu[sizeb]/I

ntrace2simu[sizeb]/I

ntrace3simu[sizeb]/I

Ntrace4simu[sizeb]/I

chi2btkb[sizeb]/F frank.meisel@lhep.unibe.ch

Numbereventb[]

showerID/I

isizeb/l

ntrace1simu[sizeb]/I

ntrace2simu[sizeb]/I

ntrace3simu[sizeb]/I

Ntrace4simu[sizeb]/I

chi2btkb[sizeb]/F

11

Reconstruction: Mode A

ShowerInstance = new EdbShowerRec();

 Choose Initiator Basetracks:

ShowerInstance-> SetInBTArray(TObjArray* InBTArray);

- Set the PVR object on ShowerInstance-> SetEdbPVRec(EdbPVRec* Ali); which the reconstruction shall operate:
- Start Reconstruction

ShowerInstance-> Execute();

Retrieve Reconstructed
 ShowerInstance-> GetRecoShowerArray();

Reconstruction: Mode B Frederics Algo Implementation

- To be filled and explained...:
- void recdown(int num, int MAXPLATE, int DATA, int Ncand, double* x0, double* y0, double* z0, double* tx0, double* ty0, double* chi20, int* W0, double* P0, int* Flag0, int* plate0, int* id0, int* TRid, double* Esim, int piece2, int piece2par)

Shower Reconstruction Algorithms

- There have been tested several algorithms to find out optimal shower reconstruction.
- Some are similar to each other, some are different.
- They can be grouped into two types:
 - Local Basetrack Cut Algorithms (CT, CA, OI, RC, BW)
 - Global Basetrack Cut Algorithms (CL, TC, GS)
- Pseudocodes of all algorithms can be found in thesis, code for OI and GS algo are given in next two slides.
- Default Used Algorithms will be Algo

OI and GS

Algorithm Shower reconstruction: CT (ConeTube) Require: BT_{IN} . Volume of plates containing basetracks. Ensure: PARA0: ConeAngle, PARA1: TubeRadius. Ensure: PARA2: Connection ΔR , PARA3: Connection $\Delta \hat{\Theta}$. Create small volume around BT_{IN} with $\pm \Delta X = \pm \Delta Y = 1200 \mu m$ $plate = plate(BT_{IN})$ while plate!= last plate do for i = 0 to $N_{BT}(Pl(i))$ do Check BT for criteria: ConeTube: if $\ni ConeTube(PARA0.PARA1) BT(i, IN)$ then CONTINUE end if for j = 0 to $N_{BT}(Shower)$) do Check BT for criteria: Connection: if $plate(BT_{IN}) != plate(BT_i)$ then if $\Delta R BT(i,j) > PARA2$ or $\Delta \Theta BT(i,j) > PARA3$ or $\Delta plate BT(i,j) > 3$ then CONTINUE end if end if if $plate(BT_{IN}) == plate(BT_i)$ then if $\Delta R BT(i, IN) > 50 \mu m$ or $\Delta \hat{\Theta} BT(i, IN) > 80 mrad$ then CONTINUE end if end if end for end for Add BT(i) to Shower Next (Ascending Z) Plate end while

return Collection of BTs = Shower

```
Algorithm Shower reconstruction: GS (GammaSearch) — !!! PSEUDOCODE !!! —
Require: BT_{IN}, Vertex VTX, Volume of plates containing basetracks.
Ensure: PARA0: Connection \Delta R_{IP}, PARA1: Connection \Delta R_{minDist}, PARA2: Connection \Delta R_{spatial},
Ensure: PARA3: Connection \Delta Z, PARA4: \Delta \Theta, PARA5: Plate Difference \Delta N_{pl},
  Create small volume around BT_{IN} with \pm \Delta X = \pm \Delta Y = 1200 \mu m
  plate = plate(BT_{IN})
  while plate!= last plate do
     for i = 0 to N_{BT}(Pl(i)) do
        Check BT for criteria: IP Distance:
        if Min(IP(BT(i), VTX), IP(BT(IN), VTX)) > PARAO then
          CONTINUE
        end if
        Check BT for criteria: Distance:
        if \Delta R_{i,IN;backward}(BT(i,IN)) > PARA2 then
          CONTINUE
        end if
        Check BT for criteria: Minimum Distance:
        if \Delta R_{minDist}(BT(i, IN)) > PARA1 then
          CONTINUE
        end if
        Check BT for criteria: Angular Distance:
        if \Delta\Theta(BT(i,IN)) > PARA4 then
          CONTINUE
        end if
        Check BT for criteria: Z Distance:
        if \Delta Z(BT(i, VTX)) > PARA3 then
          CONTINUE
        end if
        Check BT for criteria: Plate Difference:
        if Abs(plate(BT_{IN} - plate(BT_i))) > PARA5 then
          CONTINUE
        end if
     end for
     Add BT(i) to Shower
     Next (Ascending Z) Plate
  end while
  return Collection of BTs = Shower
```

- Tries to find Gamma Basetrack Pairs, i.e. One basetrack coming from e- and one from e+
- Standard Scanning way is to eyecheck a double BT pair and to identify it on the screen.
- What if we wanna search full volume? → Need automatisation!
- Usually they are in the same emulsion, but inefficiency (if not manually checked) has to be taken into account
- This is what we are trying to do here.

Two important scenarios have been identified:

SCENARIO A:

 Vertex is already located. Gammas coming from the event should point to this vertex.

•

SCENARIO B:

 Unknown vertex/event topology. Do a "blind" search over the full volume for BT pairs.

- Two important scenarios have been identified:
- SCENARIO A:

true vertex

delta Z variable

• SCENARIO B:



helper vertex

delta Z fix

- Two important scenarios have been identified:
- SCENARIO A:
 - Vertex is already located. Gammas coming from the event should point to this vertex.
 - EASY WAY, FAST, EXPECTED BG LOW
- SCENARIO B:
 - Unknown vertex/event topology. Do a "blind" search over the full volume for BT pairs.
 - DIFFICULT WAY, VERY SLOW, EXPECTED BG HIGH

- Two important scenarios have been identified:
- SCENARIO A: ← PREFERRED METHOD
 - Vertex is already located. Gammas coming from the event should point to this vertex.
 - EASY WAY, FAST, EXPECTED BG LOW
- SCENARIO B:
 - Unknown vertex/event topology. Do a "blind" search over the full volume for BT pairs.
 - DIFFICULT WAY, VERY SLOW, EXPECTED BG HIGH

Results for Shower Reco Algs:

• OI • GS

Shower Energy Measurement

See extra file...

Shower ID

Problem:

 One single Shower can originate from either photons, electrons or (charged) pions

Solution:

Shower ID class, which does a differentiation!

Technicalities:


The ID class determines Shower from its originating particle

Shower ID


- Instead of assigning the ANN one number for
 - Gamma, Electron, Pion(+-)
- We do decions based on "X vs Y", i.e.
 - Gamma vs Electron
 - Electron vs Pion
- Because the ANN is better in discriminating two-by-two variable discrimination.
- In continuation with former implementations, we save these results in eProb1,eProb90_X_vs_Y

Shower ID

- Two type of variables:
 - SHID: 0(G),1(E),2(P)
 - SPID: 0,1,2
 - So we can assign a unique
 - ID for each combination:
 - Example:
 - gamma vs e: SHID=0,SPID=1
 - e vs pi: SHID=1,SPID=2
 - 0 if it belongs to Type 0 (SHID) Always put
 - 1 if it belongs to Type 1 (SPID) frank.meisel@lhep.unibe.ch Always put



Block diagram: ID



More here...

The Cleaning of High BG level...


- We try to keep a constant BT density level of not more than 40BT/mm2
- To get this, we have the EdbPVRec patterns object, which we use for checks:
- We calculate BT density plate by plate and can adapt the Quality Cut plate by plate accordingly.


EdbPVRQuality


- This class does BG cleaning... How it works:
- EdbPVRQuality* Qualityobject = new EdbPVRQuality(gAli);
 - gAli is the filled EdbPVRec object. Usually this implicates the data coming from cp.root files.
 - LinkedTracks.root EdbPVRec object usually is already very "clean", i.e. that linked tracks are often "good" (physically) objects.
- Qualityobject->Print();
- Qualityobject->Help();

EdbPVRQuality Remove FakeDoublets.


- From Scanning we know that by view overlap of the microscope there is fake double Bts:
 - One BT scanned in two views appears (due to Sysal bending correction function) as two slightly shifted different Bts!
 - This we dont want, since several reasons. One of them is for the GammaPairSearch: The pair algo might take these as a e+e—Pair.
- Those duplicated tracks show a distinct behaviour in the dR-dT plane.
- Plot: photon pairs and instrumental BG doublets in ther dR-dT plane
- EdbPVRQuality* Qualityobject = new EdbPVRQuality(gAli);
- Qualityobject->Remove_DoubleBT();
- Returns a new EdbPVRec Object, where these doublets are eliminated.









START RECO







How the background cleaning is done.
Default Configuration.

Cleaned Volume (EdbPVRec)